



CRO

B cell Select[®] Overview – Anti-Idiotypic Antibodies



IpA


IMMUNOPRECISE ANTIBODIES

ENGINEERED for the scientific RACE

End-to-(no) end

Our mission is to provide a HUB of the most advanced intelligence and technology to treat disease, bar none.

Our goal is to improve the specificity of biotherapeutics by unlocking the language of the genome.



Anti-Idiotypic Antibodies

Superiority of Rabbit B Cell Select[®]

IPA's proprietary Rabbit B cell Select[®] workflow provides an anti-idiotypic solution that is unmatched in the industry, with a 100% program success rate in generating idiotype specific antibodies.



Anti-Idiotype Antibodies


Superiority of Rabbit B Cell Select[®]

IPA's ability to combine the high specificity and affinity of the rabbit's immune system with the unbiased depth of our "Function-First" B cell Select[®] workflow, results in the early identification of a superior panel of anti-idiotype antibodies.

25

Recent or in progress
Rabbit B cell Select[®] Anti-ID Programs

Immunogens



100%¹

Success of completed
programs from a variety of
immunogen formats

1. See data slide #6

Anti-ID

Typical program objectives



1. Program Goals

- ✓ To generate anti-idiotypic antibodies with specificity to a number of antibody formats including:
 - Fab
 - Whole IgG
 - scFv
 - VHH
 - CAR's
 - ADC's

2. Desired Properties

- ✓ Must be able to discriminate target mAb within a matrix of human immunoglobulin (human serum)
- ✓ Must not bind outside the idiotype of the antibody or cross react with other antibodies with same scaffold
- ✓ High affinity
- ✓ Both blocking and non-blocking properties
- ✓ Antibody pairs


3. Ab Development

- ✓ **B cell Select**[®]
 - Immunization: 2-5 x rabbit
 - Functional Screening: ELISA, Octet[®]
 - Timeline: ~2.5 months
- ✓ Recombinant Production

Anti-ID

Rabbit anti-ID program and screening workflow

- Pre-Phase** ● **Target Validation:** QC of immunogen prior to immunization
- Phase I** ● **Immunization:** 28-day or 78-day immunization of cohort of 2 x rabbits
- Phase II** ● **PBMC isolation:** 30 mL whole blood (no sacrifice required)
● **B cell enrichment:** Proprietary antigen specific enrichment
● **B cell culture and primary screen:** 40x 96-well plates
- Phase III** ● **Molecular cloning and sequencing:** up to 48 lead candidates based on primary screening data
● **Secondary screening:** Sequence confirmation and further testing of up to 48 recombinant lead candidates
- Phase IV** ● **Recombinant expression and purification:** up to 48 clones at scale desired by client
● **Tertiary testing:** Testing of purified material



Target Validation

- QC by ELISA, Octet®



Immunization:

- 28 day
- 78 day





Test Bleed:

- ELISA: Immunogen
- ELISA: Off-target



1° Screening (B cell supernatant):


- ELISA: Target
- ELISA: Off-target scaffold
- ELISA: Target human serum



Octet®:


- Binding
- competition
- Off-rate

IPA can determine blocking/nonblocking and relative off-rates of lead candidates prior to sequencing.



2° / 3° Screening (Recombinant)

- ELISA: Target
- ELISA: Off-target scaffold
- ELISA: Target human serum
- Octet®: Binding/competition/kinetics



Anti-ID


Octet[®] characterization of primary B cell supernatant



IPA's "Function-First" B cell Select[®] workflow allows for the identification of blocking/non-blocking functionality and in addition to off-rate ranking [prior to cloning/sequencing](#)


Competition Assay

Determine blocking / non-blocking properties



Off-rate Ranking / Full Kinetics

Relative & absolute affinity measurements



Anti-ID


Completed anti-Id rabbit programs since 2020



B cell Select® Rabbit Anti-ID Programs

Program Average sorted by		Normalized to 3840 clone screen		Normalized to cloning and sequencing of top 48 clones		
		Total Ag (+) Clones <i>Primary Screen</i>	ID-Specific Clones <i>Secondary Screen</i>	B Cell Cloning Efficiency	Unique Sequences	
					<i>Combined H + L</i>	<i>CDR3</i>
Immunogen	Fab	160 5%	91 57%	37 77%	36 97%	33 88%
	scFv	101 3%	57 56%	38 80%	38 99%	36 94%
Immunization	Rapid Prime (28-Day)	147 4%	72 49%	43 80%	42 97%	35 83%
	Standard (78-days)	135 4%	85 63%	34 69%	34 99%	33 98%
Consolidated Averages		142 4%	67 47%	35 72%	34 97%	31 89%


Immunization Method Breakdown



100%

Program success in generating
idiotype specific antibodies

Immunogen Breakdown



B cell Select®

Consolidated program data average from last 2 years




All B cell Select® Programs since 2020 including Anti-ID of all Species

B cell programs	Total Ag (+) Clones <i>Primary Screen</i>	B Cell Cloning Efficiency	Unique Sequences	
			<i>Combined H + L</i>	<i>CDR3</i>
66	7%	82%	90%	62%
Completed or in Progress	ELISA & Flow	Transgenic, Rabbit, Canine, Chicken, & Rodent	One (1) A.A. difference in stitched heavy and light chain within a program	One (1) A.A. difference in CDR3 within a program


Rodent Strain Breakdown



Program Breakdown



Target Breakdown





Engineered for the race
and the shared pursuit
of clinical success



CONTACT

Dr. Barry Duplantis
VP of Client Relations

bduplantis@immunoprecise.com

www.immunoprecise.com

IPA-Canada

Vancouver Island Technology Park
Unit 3204-4464
Markham Street
Victoria, British Columbia,
V8Z 7X8, Canada

IPA-Oss

Pivot Park RE2142
Industrielaan 63
5349 AE Oss
The Netherlands

IPA-Utrecht

Utrecht Science Park
Life Science Incubator
Yalelaan '62
3584 CM Utrecht
The Netherlands



IMMUNOPRECISE ANTIBODIES