


Highly diverse human scFv libraries


The Challenge

Neoantigens: promising, but challenging targets for developing cancer therapeutics. Due to their tumor cell-restricted expression, neoantigens are attractive targets for cancer therapy. Most neoantigens derive from intracellular proteins, making the intact neoantigen a difficult target for antibody-mediated therapy. However, neoantigens, like most endogenous proteins, are processed by the proteasome, and resulting neoantigen peptides may form complexes with HLA type 1 and 2 proteins, called neopeptides, that locate to the plasma membrane for presentation to T-cells (figure on right) and are therefore potentially accessible to antibodies. However, due to the minor differences between neopeptides and their non-mutated counterpart it is challenging to obtain neopeptide-specific antibodies through immunizations with recombinant neopeptides. Rather, it requires antibody discovery techniques that allow for stringent selection procedures, such as phage display. IPA's highly diverse human scFv phage libraries have proven to be successful in obtaining antibodies against difficult targets and have the advantage of containing antibodies with favorable developability profiles. IPA's human phage libraries could therefore be an excellent source for clinically-developable neopeptide-specific antibodies. Here, we present a case study identifying antibody fragments (scFvs) against a common neopeptide by utilizing IPA's in-house human scFv libraries.




The Solution

Tackling neoantigen antibody discovery challenges with fully human scFv libraries




IPA-EU's human scFv phage libraries with favorable developability characteristics (figure on left) were subjected to five consecutive panning rounds consisting of negative selections against wildtype and irrelevant peptide HLA (pHLA) complexes and a positive selection against the neoepitope pHLA complex (figure on right). Selection and washing stringencies were increased after each round.

Output summary: accelerated neopeptide antibody discovery




Polyclonal panning outputs from each selection round were screened for reactivity towards the neopeptide, wildtype and irrelevant pHLA complex and BSA as negative control in an ELISA. Neopeptide differentiating reactivity is apparent from selection round 3 on and increases until selection round 5. Round 4 and 5 panning outputs were selected for monoclonal scFv analysis.



ELISA-based reactivity analysis (figure above) and subsequent sequencing of monoclonal scFvs yielded clones with a range of specificities (from neopeptide specific to cross-reactive with all three pHLA complexes) that belonged to 12 unique sequence families, labeled A-L. Four of these families showed highly specific neopeptide pHLA complex reactivity. Interestingly, these four families (A, B, C and I) cluster closely as shown in a phylogenetic tree analyzing full VH and VL sequences (figure on right).

Streamlined workflow for rapid confirmation of neopeptide specificity



Each unique scFv was recombinantly expressed and purified, followed by screening for reactivity towards neopeptide, wildtype and irrelevant pHLA complex and BSA as negative control in an ELISA. In general, reactivity profiles were confirmed with some deviations (Family 2 and Family 11).

Conclusion

- IPA's highly diverse human scFv libraries were successfully used to identify neopeptide targeting antibodies, yielding 11 VH-CDR3 families, of which four showed highly neopeptide specific reactivity.
- Recombinant production and screening of scFvs confirmed the identification of 4 neopeptide specific antibodies.
- Phage display-mediated antibody discovery using fully human repertoires allows for rapid and accurate selection of specific, clinically-suitable neopeptide antibodies.

For more information — email: info@ipatherapeutics.com web: ipatherapeutics.com



UNITED STATES | CANADA | EUROPE

©2023 ImmunoPrecise Antibodies LTD. All rights reserved. IPA is a trademark of ImmunoPrecise Antibodies LTD.